



P H Y S I Q

THREE SIXTY°

**A REVOLUTION IN
BODY CONTOURING**

How Do Non-Invasive Body Treatments Work?



**TEMPERATURE
CONTROLLED**



**MUSCLE
STIMULATION**



**MULTI-
MODALITY**



What is PHYSIQ 360?

- PHYSIQ 360 uses 4 applicators to target stubborn areas where diet & exercise may not be enough
- 3 energy modalities:
 - Burn fat with laser energy (LZR)
 - Stimulate muscles with EMS
 - Do both in one session with STEP (Sequential Thermal & Electrical Pulse) Technology
- Fast and comfortable treatments
- No downtime

Treat Fat & Muscle



HEAT DELIVERED TO FAT

PURE BEAM Technology

Delivers heat to deeper layers of
fat

Heats a larger volume of tissue
Maximize comfort



COMPLETE MUSCLE RE-EDUCATION

Electrical Muscle Stimulation

Re-educate muscle fibers
Engages muscles in ways
working out can't




STEP Technology



UNIQUE STEP TECHNOLOGY

Alternates between
burning fat & stimulating
muscles in one session

Maximizes results without
increasing treatment time



Total Body – *Every* Body



PHYSIQ 360 is the most complete body contouring solution on the market:


- Target fat and muscle
- Applicators fit large & small body areas
- Treat multiple body areas at once
- Safe & effective for all skin/body types
- Complements medical weight loss
- Built-in lymphatic drainage PLUS handpiece

What is Lymphatic Drainage Massage (LDM)?



- LDM is a proven method of increasing circulation and can help expedite the natural healing process
- PHYSIQ 360 kills the fat cells – LDM helps your lymphatic system process the fat cells out of the body
- PHYSIQ 360 has a special handpiece for LDM treatments
- Treatments are fast (2 minutes/area) and completely comfortable

PHYSIQ 360 Treatment Areas



PHYSIQ 360 Treatment – What to Expect



Up to 4 applicators will be placed on treatment area(s)




Treatment will last 15-36 minutes



EMS feels like waves of contractions or sit-up/crunch



LZR feels like a warm stone massage



A 2-minute lymphatic drainage massage is performed post-treatment



Treatments are performed every 7-14 days



Average of 5 treatments are recommended



Pre-Treatment



Avoid direct sun exposure (2-4 weeks)



Suspend Accutane (4-6 months)



Proper hydration (3-4 days)



Treatment area clean shaven



Comprehensive consultation

Post-Treatment



Self-massage 2x daily with PHYSIQ lotion



Increase water intake



Decrease alcohol consumption



Gentle exercise is recommended



Avoid direct sun exposure



Hear from Our Patients

“I’m in and out within 45 minutes.”

“Honestly, I fall asleep every time!”

“I LOVE the feeling of the muscle contractions; I feel like I’m getting a workout laying down!”

“I feel so much more confident in my clothes.”



“Helped me with that one stubborn area that wouldn’t budge with anything else.”

THREE SIXTY°